

Technical Report 5-32902  
Contract No. NAS8-38609  
Delivery Order No. 31  
UAH Research Report No.

*NAS8-38609 D.O. 31*

Vacuum Chamber and Adaptors for  
Shearography and Holography  
(5-32902)

Final Technical Report for Period  
8 April 1992 through 8 January 1993

January 1993


Prepared by

James R. Clark  
Charles M. Horton  
Robert J. Campbell

Research Institute  
The University of Alabama in Huntsville  
Huntsville, Alabama 35899

Prepared for

George C. Marshall Space Flight Center  
National Aeronautics and Space Administration  
Marshall Space Flight Center, Alabama 35812  
Attn: EH-13 (Dr. Samuel S. Russell)

NASA National Aeronautical and Space Agency		Report Document Page	
1. Report No.  Final Technical	2. Government Accession No.	3. Recipient's Catalog No.	
4. Title and Subtitle  "Vacuum Chamber & Adaptors for Shearography & Holography"		5. Report Due  January 1993	
		6. Performing Organization Code  Research Institute, UAH	
7. Author(s)  James Clark Charles Horton		8. Performing Organization Report No.	
9. Performing Organization Name and Address  UAH Research Institute RI E-47 Huntsville, AL 35899		10. Work Unit No.  5-32902	
		11. Contract or Grant No.  NAS8-38609 D.O. 31	
12. Sponsoring Agency Name and Address  National Aeronautics and Space Administration Washington, D.C. 20546-001 Marshall Space Flight Center, AL 35812		13. Type of report and Period covered  Final Technical Report	
		14. Sponsoring Agency Code	
15. Supplementary Notes			
16. Abstract  Study methods to excite aerospace components using vacuum to induce deformations which will be measured by shearography, holography, and image correlation.			
17. Key Words (Suggested by Author(s))  Vacuum Excitation Chamber		18. Distribution Statement  <i>Public Avail per telecon</i>  <i>w/ Joyce Turner</i> <i>cjb 3-31-94</i>	
19. Security Class. (of this report)  Unclassified	20. Security Class. (of this page)  Unclassified	21. No. of pages  6 + Appendix	22. Price

PREFACE

This technical report was prepared by the staff of the Research Institute, The University of Alabama in Huntsville. The purpose of this report is to provide documentation of the work performed and the results obtained under the Delivery Order 31 of Marshall Space Flight Center (MSFC) Contract No. NAS8-38609. Mr. James R. Clark was the Principal Investigator and Mr. Charles M. Horton was the Co-Principal Investigator for this nine month level of effort. Dr. Samuel S. Russell of the NDE Branch, EH-13, of MSFC/NASA provided technical support and direction. Mr. Thomas J. King of Pratt & Whitney provided support in the mechanics of Holography/Shearography Systems.

The views, opinions, and/or findings contained in this report are those of the author(s) and should not be construed as an official NASA position, policy, or decision unless so designated by other official documentation.

I have reviewed this report, dated 25 Jan. 93 and the report contains no classified information.

  
Principal Investigator

  
Co-Principal Investigator

Approval:

  
Research Institute

## TABLE OF CONTENTS

	Page
1.0 INTRODUCTION .....	1
2.0 ACTIVITIES .....	1
3.0 CONCEPTUAL DESIGN .....	2
4.0 PROTOTYPE DEVELOPMENT & OPERATION .....	2
5.0 CONCLUSIONS & RECOMMENDATIONS .....	6

### APPENDICES

A. UNBONDED AREAS DETECTED BY SHEAROGRAPHY USING PRESSURE  
REDUCTION

B. VACUUM CHAMBER DESIGN DRAWINGS

C. VACUUM CHAMBER STRESS ANALYSIS

D. VACUUM CHAMBER PARTS LIST

## 1.0 INTRODUCTION

Electronic Holographic and Shearographic NDE are optical inspection methods that can detect sub-surface flaws such as unbonds and delaminations in bonded, laminated, composite or coated structures. The two methods detect flaws by displaying, in real-time on a video monitor, the surface deformation (strain) that results from the structure being stressed by either a steady state load (pressure, vacuum, or mechanical) or a dynamic load (temperature or acoustic vibration). Since the structure or coating is locally weak where a sub-surface flaw exists, Holographic or Shearographic NDE will detect an area of high surface deformation (strain) when a load is applied.

## 2.0 ACTIVITIES

NASA and UAH collaborated in a joint study to consider methods of vacuum exciting structures to induce deformations to be measured by shearography, holography, and image correlation. A Pratt & Whitney (P&W) Electronic Holography/Shearography Inspection System (EH/SIS) was used by the NDE Branch, EH-13, of NASA/MSFC in this study. The system uses a ND:YAG laser and a frequency doubler to produce a constant wave 532nm green light source for illumination of the test object. Light is reflected from the object and acquired by a standard CCD Video Camera, equipped with special optics, for image processing. Images may be acquired through any optically clear non-distorting medium.

### 3.0 CONCEPTUAL DESIGN

The Bell Jar is the most commonly recognized device for a controlled vacuum environment in scientific experiments. A cylindrical glass vessel would be the ideal choice as a vacuum chamber because the test object could be viewed from any direction without opening the chamber to rotate the piece. Due to the limitations of the Bell Jars size and the optical distortion caused by the curvature of the vessel walls another design was required.

### 4.0 PROTOTYPE DEVELOPMENT AND OPERATION

The design chosen was a hexagon shaped vacuum chamber made of six flat acrylic panels. The flat optically clear panels would provide distortion free viewing from all sides. The vacuum excitation chamber consists of three major components; a lid, chamber body, and base plate. (Figure 1.) The vacuum chamber design was scaled and constructed to hold larger size specimens than available vacuum chambers. The chamber accommodates a specimen of 1 meter in diameter (max). The chamber body was constructed of 22" x 48" x 1" acrylic panels solvent welded to form one continuous hexagon tube, open at both ends. The top and bottom edges are fitted with molded continuous neoprene rubber channel gaskets to insure a vacuum seal. The removable lid consisted of a flat hexagon of 1 1/2" acrylic material. The oversized lid rests on the channel gasket that covers the top of




Figure 1.  
Lid, Chamber Body, and Base Plate.

the vacuum chamber, closing the open end. The channel gasket that covers the bottom of the chamber rests upon a reinforced 48" x 48" x 1 1/2" aluminum base plate. The chamber has four handles so that it may be easily moved from the isolation table to its storage cart. Care should be taken not to roll the channel gasket off the bottom edge of the chamber during movement. The chamber should be lifted clear of

the base plate and not dragged when moved.

Because vacuum deformation of the base plate would cause significant movement of the test specimen, a false bottom was constructed of 1/2" plate aluminum to fit inside the test chamber. This hexagon shaped "table" fits easily inside the chamber with support legs at the perimeter. The flat plate of the table was drilled and tapped with 1/4"-20 SAE holes to allow for the installation of support members and fixtures for specimens under inspection. Any deflection that occurs during extreme vacuum conditions, is most significant in the middle of the base

plate and less evident at the chamber walls. Because the table is supported along the perimeter of the chamber walls the "free floating" table dampens any movement due to vacuum deformation.

With the lid removed the chamber could be tilted on its edge for easy insertion and removal of the specimen table under the raised chamber edge. Easy access could also be gained for any small test objects that required inspected. (Figure 2.)




Figure 2.

An Industrial wet/dry vacuum head installed on a 55 gallon drum was used to provide a vacuum source to quickly reduce the pressure in the chamber. The vacuums 1 1/2" hose pipe is attached to the vacuum out gas ball cock by a quick release connector. The 3.5 HP two-stage bypass motor moves 110 CFM of air at the static pressure rate of 7.7in HG (3.8 PSIG). This volume enabled the chamber pressure to be reduced and held at 7.7in HG (3.8 PSIG) in 18 seconds. A convectional vacuum pump requires at least 10 times longer to achieve the same vacuum. Rapid changes in the internal pressures of the test specimen with respect to the chamber environment provides more possibility of movement in unbonded materials. Slower vacuum rates allow for the equalization of the internal pressure differentials which cause detectable movement. (Figure 3.)



Once the desired pressure had been reached the vacuum out gas ball cock is closed to maintain the desired pressure in the test chamber. A second gas ball cock is installed as a "waste gate". The waste gate is adjusted until the air flow entering the chamber equaled the air flow

exiting the chamber. This "fine tuning" of the waste gate aids in establishing an equilibrium between atmosphere and the vacuum source when constant vacuum levels are required. The waste gate also allows pulsing of the vacuum level in the chamber. A needle valve is provided for precision control of air flows entering the chamber. In cases where heavier steady state loads are needed to stress a part, a 1/4" hose barb connector attached to the needle valve allows any conventional vacuum pump to be attached and further reduce the pressure in the chamber.

As a safety precaution a vacuum breaker has been installed on the chamber. The




Figure 3.  
Inlet/Outlet Valves




Figure 4.  
Vacuum Breaker

vacuum breaker has been set at 13in. HG (6.4 PSIG). Once the vacuum inside the chamber reaches the preset level the breaker opens to the outside atmosphere and stabilizes the pressure in the chamber. (Figure 4.)

## 5.0 CONCLUSIONS AND RECOMMENDATIONS


The vacuum excitation chamber has proven to be an excellent method for stressing aerospace components. The results achieved so far indicate that Holographic/Shearographic Optical NDE is a viable method for detecting surface deformations resulting from sub-surface flaws. Further evaluation would appear to be in order to validate this method as an acceptable means of NDE inspection.

EXHIBIT A


UNBONDED AREAS DETECTED BY SHEAROGRAPHY  
USING PRESSURE REDUCTION

ORIGINAL PAGE  
COLOR PHOTOGRAPH

Typical Shearographic Image



Pressure Reduction




Pressure Reduction

ORIGINAL PAGE IS  
OF POOR QUALITY


**EXHIBIT B**

**VACUUM CHAMBER DESIGN DRAWINGS**





VACUUM EXCITATION CHAMBER  
LID, CHAMBER BODY, AND BASE PLATE

COVER MADE OF 1 1/2" TK. MATERIAL




VACUUM EXCITATION CHAMBER LID






VACUUM EXCITATION CHAMBER FOOTPRINT  
TOP VIEW



CHAMBER WITH HANDELS AND HOLES



VACUUM EXCITATION CHAMBER PANEL LOCATIONS  
TOP VIEW




NOTE: HANDLES WILL BE GLUED AND BOLTED.

1/4" - 20 X 3" BOLT OR FABRICATORS CHOICE.

TAPPED HOLE MUST NOT BREAK THROUGH TO INSIDE OF CHAMBER WALL.




## VACUUM EXCITATION CHAMBER HANDLE CONSTRUCTION








ENLARGED TOP VIEW

# VACUUM EXCITATION CHAMBER PANEL 2






ENLARGED TOP VIEW

# VACUUM EXCITATION CHAMBER PANEL 3



VACUUM EXCITATION CHAMBER  
 PANEL 4






ENLARGED TOP VIEW


# VACUUM EXCITATION CHAMBER PANEL 6






MATERIALS: 0.5" ALUM. PLATE

VACUUM EXCITATION CHAMBER  
FLOATING BOTTOM




VACUUM EXCITATION CHAMBER  
FLOATING BOTTOM LEGS




### CHAMBER FOOTPRINT

11 FEET IN CIRCUMFERENCE



### ENLARGED PANEL SECTION



### ENLARGED CROSS SECTION OF CHANNEL GASKET

TOP CHANNEL GASKET 15 +/- DURO

BOTTOM CHANNEL GASKET 30 +/- DURO

GASKETS ARE MADE OF NEOPRENE RUBBER

## VACUUM EXCITATION CHAMBER CHANNEL GASKETS

EXHIBIT C


VACUUM CHAMBER STRESS ANALYSIS

# PROPERTIES OF CAST ACRYLIC SHEET

.250 INCHES THICK - ASTM TEST METHOD

SPECIFIC GRAVITY (D792) .....	1.19
TENSIL STRENGTH (D638) .....	10,000psi (69MPA AT 4.2%)
ELASTIC MODULUS .....	400,000psi (2800 MDA)
FLEXURAL STRENGTH (D790) .....	16,500psi (RUPTURE 114 MPA)
MODULUS OF ELASTICITY .....	475,000psi (3300 MPA)
COMPRESSION STRENGTH (D695) .....	18,000psi (YIELD 124 MPA)
MODULUS OF ELASTICITY .....	430,000psi
SHEAR STRENGTH (D732) .....	9,000psi
IMPACT STRENGTH (D256) .....	0.4 lbs./in. of notch (IZOD MILLED NOTCH) (21.6j/m of notch)
ROCKWELL HARDNESS (D785) .....	M94
BARCOL HARDNESS (D2583) .....	49
RESIDUAL SHRINKAGE .....	2% per 0.250 inch (INTERNAL STRAIN)
OPTICAL CLEAR MATERIAL	
REFRACTIVE INDEX (D542) .....	.025 - 1.49
LIGHT TRANSMISSION (D1003) .....	92%
UV TRANSMISSION .....	0 at 320nm
HAZE .....	LESS THAN 1%
THERMAL .....	300deg.F - 350deg.F
WATER ABSORPTION .....	24 hours - 73deg.F - 2%

Stress loads for a single 22" x 48" x 1" vacuum chamber wall plate. Using the stress formulas for flat plates with straight



boundaries and constant thickness the maximum safe load limits have been calculated. The tensile strength of acrylic is 10,000psi. With a safety factor of five the maximum safe psi load levels for a flat acrylic plate should not exceed 2,000psi.

$$\sigma = \frac{\beta \times b^2 \times q}{t^2}$$


$$a = 48'' \quad b = 22'' \quad \frac{a}{b} = 2.18 \quad \beta = 0.6102 \quad q = \text{lbs per inch}^2$$

$$\sigma = \frac{0.6102 \times (22 \text{ inch})^2}{(1 \text{ inch})^2} = q \times 295.3 \text{ psia}$$

$$\sigma = q \times 295.3 \text{ psia}$$

$$\text{For } q = 6.5 \text{ psia} \quad \sigma = 1919.45 \text{ psia}$$

Stress loads for the 38" x 44" x 1 1/2" vacuum chamber lid. Using the stress formulas for flat plates with straight



boundries and constant thickness, the maximum safe load limits have been calculated. The tensil strength of acrylic is 10,000psi. With a safety factor of five the maximum safe psi load levels for a flat acrylic plate should not exceed 2,000psi.

$$\sigma = \frac{\beta \times b^2 \times q}{t^2}$$


$$a = 44'' \quad b = 38'' \quad \frac{a}{b} = 1.157 \quad \beta = 0.3762 \quad q = \text{lbs per inch}^2$$

$$\sigma = \frac{0.3762 \times (38 \text{ inch})^2}{(1.5 \text{ inch})^2} = q \times 241.43 \text{ psia}$$

$$\sigma = q \times 241.43 \text{ psia}$$

$$\text{For } q = 6.5 \text{ psia} \quad \sigma = 1569.33 \text{ psia}$$

Approximation of stresses for vacuum chamber using the



formulas for membrane stresses and deformations in thin-walled vessels. These number represent a cylinder 44" in diameter with 1" thick walls. The tensile strength of acrylic is 10,000psi. With a safety factor of five the maximum safe load limits have been calculated and should not exceed 2,000psi.

$$\sigma = \frac{-q \times r}{t}$$

$$t = \text{wall thickness} = 1''$$

$$q = \text{lbs/inch}^2$$

$$r = \text{radius} = 22''$$


$$\sigma = \frac{-1 \times 22''}{t}$$

$$\sigma = 1 \text{ lb} \times -22'' = -22 \text{ psia}$$

$$\sigma = 6.5 \text{ lb} \times -22'' = -143 \text{ psia}$$



Approximation of stresses for vacuum chamber using the



formulas for membrane stresses and deformations in thin-walled vessels.

These numbers represent a cylinder 44" in diameter with 1" walls. The tensile strength of acrylic is 10,000psi. With a safety factor of five the maximum safe load limits have been calculated and should not exceed 2,000psi.

$$\sigma = \frac{P}{t}$$

$t = \text{wall thickness} = 1''$

$p = \text{lbs/linear inch}$

$$cir = 2 \times \pi \times radius = 2 \times 3.14 \times 22'' = 138.23''$$

$$\sigma = 1lb \times 138.23'' = 138.23psia$$

$$\sigma = 6.5lb \times 138.23'' = 898.49psia$$

EXHIBIT D

VACUUM CHAMBER PARTS LIST

## PARTS LIST - VACUUM EXCITATION CHAMBER

1. VACUUM EXCITATION CHAMBER - 1" ACRYLIC PLATE  
(CHAMBER WEIGHT - 250 LBS.)
2. CHAMBER LID - 1 1/2" ACRYLIC PLATE  
(LID WEIGHT - 90 LBS.)
3. TOP EDGE GASKET ( NEOPRENE RUBBER 15 +/- 5 DURO )
4. BOTTOM EDGE GASKET ( NEOPRENE RUBBER 30 +/- 5 DURO )
5. VACUUM GAUGE, BOURDON TUBE TYPE, 4 1/2" DIAL, 1/4" PIPE SIZE
6. NEEDLE VALVE, 1/4" PIPE SIZE
7. 1/4" CLOSE NIPPLE, BRASS
8. 1/4" x 1/4" HOSE BARB, BRASS
9. 1/2" - 1/4" REDUCER, BRASS
10. 1/4" - 45DEG. ELBOW, BRASS
11. 1/8" MUFFLER/FILTER
12. 1/2" MUFFLER/FILTER
13. 1/4" TEE, BRASS
14. VACUUM BREAKER, ADJUSTABLE, 1/2" PIPE SIZE
15. 2 EA. 1 1/2" x 1/16" RUBBER WASHERS
16. 2 EA. 7/8" x 1/16" RUBBER WASHERS
17. 2 EA. 1 1/2" PVC ACID DRAINS
18. 2 EA. 1/2" BRASS DRAINS
19. 2 EA. 1 1/2" GAS BALL COCKS
20. INDUSTRIAL WET/DRY VACUUM HEAD ASSEMBLY
21. 55 GALLON FRH DRUM
22. DRUM DOLLY, 55 GALLON
23. 1 1/2" QUICK CONNECT COUPLER W/ HOSE CLAMPS

PARTS LIST CONT.

- 24. VACUUM CHAMBER STORAGE CART
- 25. 48" x 48" x 1 1/2" ALUMINUM BASE PLATE  
(BASE PLATE WEIGHT - 345 LBS.)
- 26. CHAMBER FALSE FLOOR, 1/2" PLATE ALUMINUM - HEXAGON SHAPE

